



# Persona RF PREMIUM+

*Radiographic/Fluoroscopy & Digital Radiography  
multipurpose room*



**Compact, high functionality dRF**  
dual-purpose room performs fluoroscopy &  
general radiography in one.  
Auto-positioning table & tube make it  
possible to utilize either in-room or the  
remote control console.

Scalability to add an overhead  
tube and wall bucky integrated  
with FDR D-EVO III, all with one  
generator. Minimizes spacing  
requirements and costs.



## Product highlights

- DR detector based system, for enhanced dose efficiency, speed and field of view
- 71" SID accommodates chest and long length stitching exams upright and supine
- Low table elevation enables more exam types with easier patient access.
- Clinician access behind table for easier patient assist vs. traditional systems
- Video camera to verify positioning
- Auto-positioning table and tube
- Solid construction 584 lbs. capacity without movement limitations
- Ceiling mounted in-room monitors (option)