



# Persona RF PREMIUM

*Radiographic/Fluoroscopy & Digital Radiography  
multipurpose room*

## **Compact, high functionality dRF**


dual-purpose room performs fluoroscopy & general radiography in one.

Auto-positioning table & tube make it possible to utilize either in-room or the remote control console.

Scalability to add an overhead tube and wall bucky integrated with FDR D-EVO II, all with one generator. Minimizes spacing requirements and costs.

## **Product highlights**

- DR detector based system, for enhanced dose efficiency, speed and field of view
- 71" SID accommodates chest and long length stitching exams upright and supine
- Low table elevation enables more exam types with easier patient access.
- Clinician access behind table for easier patient assist vs. traditional systems
- Video camera to verify positioning
- Auto-positioning table and tube
- Solid construction 584 lbs. capacity without movement limitations
- Ceiling mounted in-room monitors (option)



## **FUJIFILM Healthcare Americas Corporation**

81 Hartwell Avenue, Suite 300, Lexington, MA 02421  
fujifilmhealthcare.com 800.431.1850

© FUJIFILM Healthcare Americas Corporation DOC-0057429-A

**FUJIFILM**  
Value from Innovation